




Professional Manufacturer of Endoscopic Surgical Instruments





Company Introduction

01/02

Product Introduction

- 03
- 04 Disposable Non-Vascular Lumen Guide Wire
- 05 Disposable Biopsy Forceps
- 06 Disposable Digestive Guide Wire
- 07 Disposable Spray Catheter
- 08 Disposable Stone Retrieval Basket
- 09 Disposable Cytology Brush
- 10 Disposable Ureteral Access Sheath
- 11 Disposable Nasobiliary Drainage Tube
- 12 Disposable Urinary Stone Retrieval Basket
- 13 Disposable Cleaning Brush
- 14 Disposable Polypectomy Hot Snare
- 15 Disposable Hemostasis Clips
- 16 Disposable Injection Needle
- Disposable Polypectomy Cold Snare




Jiangxi Zuoruihua Medical Instrument Co., Ltd. is a modern enterprise dedicated to the research and development, production and sales of endoscopic minimally invasive surgical instruments.

In the face of the sustained and rapid growth of the global market of endoscopic minimally invasive diagnosis and treatment consumables and microwave ablation market, as well as the deepening of China's medical reform and the urgent need of countries in the world to reduce medical costs, we are seizing the new historical opportunity to accelerate the marketing of new strategic products such as OCT-3D, endoscopic early cancer diagnosis and treatment products, endoscopic ultrasound products and a new generation of microwave ablation devices.

We have obtained ISO13485 quality system certificate and CE certificate. We provide OEM service for customers and strive to become a famous brand in the world medical industry.



ISO13485




C E

Disposable Non-Vascular Lumen Guide Wire

Features

- Impenetrable soft head, fully developed under X-ray
- Double color alternate spiral pattern, the movement can be seen under the endoscope.
- Super smooth hydrophilic head, small friction and smooth intubation.
- Super elastic Ni-Ti alloy mandrel with excellent push and pass ability.
- Hydrophilic coating on the surface significantly reduces the friction.
- Smooth head end prevent mucosal tissue damage.




Specification

REF	Diameter	Working Length (mm)	Tip Style
Zebra Guide Wire	0.018" / 0.022" / 0.025" / ★ 0.028" / ★ 0.032" / ★ 0.035" / 0.038"	800 / 1200 / 1500	Straight/Angle
Loach Guide Wire	0.018" / 0.022" / 0.025" / ★ 0.028" / ★ 0.032" / ★ 0.035" / 0.038"	800 / 1200 / 1500	Straight/Angle

Disposable Biopsy Forceps

Features

- Complete specifications meet various clinical needs
- High occlusal force makes sampling safer and more effective.
- Needle: accurate positioning, easy to fix
- Jaws: easy to open and close, not stuck
- Spring Tube: super smooth PE coating, reduce the damage to clamp channel.
- Handle: simple and convenient operation



Specification

No.	REF	Jaw Open Size (mm)	Maximum Working Channel (mm)	Working Length (mm)	Serrated Jaw	Spike	PE Coating
1	ZRH-BFA-2423-CZS	6/5/4	2.4/1.8	2300/1800/1200/ 1000/600	Yes	Yes	Yes
2	ZRH-BFA-2423-CZL		2.4/1.8	2300/1800/1200/ 1000/600	Yes	Yes	No
3	ZRH-BFA-2423-CWS		2.4/1.8	2300/1800/1200/ 1000/600	Yes	No	Yes
4	ZRH-BFA-2423-CWL		2.4/1.8	2300/1800/1200/ 1000/600	Yes	No	No
5	ZRH-BFA-2423-PZS		2.4/1.8	2300/1800/1200/ 1000/600	No	Yes	Yes
6	ZRH-BFA-2423-PZL		2.4/1.8	2300/1800/1200/ 1000/600	No	Yes	No
7	ZRH-BFA-2423-PWS		2.4/1.8	2300/1800/1200/ 1000/600	No	No	Yes
8	ZRH-BFA-2423-PWL		2.4/1.8	2300/1800/1200/ 1000/600	No	No	No

Disposable Digestive Guide Wire

Features

- Impenetrable soft head, fully developed under X-ray
- Triple protection design of hydrophilic head end and inner core
- Zebra smooth coating has good trafficability and no irritation
- Anti-twist Ni-Ti alloy inner core provides excellent torsion and pushing force
- Super elastic Ni-Ti alloy mandrel with excellent push and pass ability
- Tapered design head enhance the flexibility intubation and operation success rate
- Smooth head end prevent mucosal tissue damage.




Specification

No.	REF	Straight/Angle	Maximum O.D.		Working Length ± 50 (mm)
			±0.004 (inch)	±0.1mm	
1	ZRH-XBM-W-2526	Angle	0.025	0.63	2600
2	ZRH-XBM-W-2545	Angle	0.025	0.63	4500
3	ZRH-XBM-Z-2526	Straight	0.025	0.63	2600
4	ZRH-XBM-Z-2545	Straight	0.025	0.63	4500
5	ZRH-XBM-W-3526	Angle	0.035	0.89	2600
6	ZRH-XBM-W-3545	★ Angle	0.035	0.89	4500
7	ZRH-XBM-Z-3526	Straight	0.035	0.89	2600
8	ZRH-XBM-Z-3545	★ Straight	0.035	0.89	4500

Disposable Spray Catheter

Features

- High cost performance
- Easy operation
- Needle Tube: large flow, fully reduce injection resistance
- Outer Sheath: smooth surface and smooth intubation
- Inner Sheath: smooth lumen and smoother liquid delivery
- Handle: portable single hand control




Specification

No.	REF (Injection Type)	Sheath O.D ± 0.1 (mm)	Working Length L ± 50 (mm)	Needle Gauge G (mm)	Needle Extension Length 1±1 (mm)
1	ZRH-PZ-2418-214	∅ 2.4	1800/2300	21(0.800-0.830)	4
2	ZRH-PZ-2418-234	∅ 2.4	1800/2300	23(0.600-0.673)	4
3	ZRH-PZ-2418-254	∅ 2.4	1800/2300	25(0.500-0.530)	4
4	ZRH-PZ-2418-216	∅ 2.4	1800/2300	21(0.800-0.830)	6
5	ZRH-PZ-2418-236	∅ 2.4	1800/2300	23(0.600-0.673)	6
6	ZRH-PZ-2418-256	∅ 2.4	1800/2300	25(0.500-0.530)	6
No.	REF (Spray Type)	Sheath O.D. (mm)	Working Length L ± 50(mm)	/	/
1	ZRH-PW-1810	∅ 1.8	1000	/	/
2	ZRH-PW-1818	∅ 1.8	1800	/	/
3	ZRH-PW-2418	∅ 2.4	1800	/	/
4	ZRH-PW-2423	∅ 2.4	2400	/	/

Disposable Stone Retrieval Basket

Features

- Easy operation
- Various models can be customized
- Injectable contrast agent
- Ni-Ti alloy material
- Conforming to the structure of biliary tract



Specification

No.	Basket Type	REF	Basket Diameter (mm)	Basket Length (mm)	Working Length (mm)	Channel Size (mm)	Contrast Agent Injection (Yes / No)
1	Diamond Type (A)	ZRH-BA-1807-15	15	30	700	Φ 1.9	No
2		ZRH-BA-1807-20	20	40	700	Φ 1.9	No
3		ZRH-BA-2416-20	20	40	1600	Φ 2.5	Yes
4		ZRH-BA-2416-30	30	60	1600	Φ 2.5	Yes
5		ZRH-BA-2419-20	20	40	1900	Φ 2.5	Yes
6		ZRH-BA-2419-30	30	60	1900	Φ 2.5	Yes
7	Oval Type (B)	ZRH-BB-1807-15	15	30	700	Φ 1.9	No
8		ZRH-BB-1807-20	20	40	700	Φ 1.9	No
9		★ ZRH-BB-2416-20	20	40	1600	Φ 2.5	Yes
10		ZRH-BB-2416-30	30	60	1600	Φ 2.5	Yes
11		ZRH-BB-2419-20	20	40	1900	Φ 2.5	Yes
12		★ ZRH-BB-2419-30	30	60	1900	Φ 2.5	Yes
13	Spiral Type (C)	ZRH-BC-1807-15	15	30	700	Φ 1.9	No
14		ZRH-BC-1807-20	20	40	700	Φ 1.9	No
15		★ ZRH-BC-2416-20	20	40	1600	Φ 2.5	Yes
16		ZRH-BC-2416-30	30	60	1600	Φ 2.5	Yes
17		★ ZRH-BC-2419-20	20	40	1900	Φ 2.5	Yes
18		ZRH-BC-2419-30	20	60	1900	Φ 2.5	Yes

Disposable Cytology Brush

Features

- Thin brushes for an optimal collection of cells
- Includes plastic tube and metal head for closure
- Sterile single packaging




Specification

No.	REF	Brush Diameter (mm)	Brush Length (mm)	Working Length (mm)	Maximum insert width (mm)
1	ZRH-CB-1812-2	Φ 2.0	10	1200	Φ 1.9
2	ZRH-CB-1812-3	Φ 3.0	10	1200	Φ 1.9
3	ZRH-CB-1816-2	Φ 2.0	10	1600	Φ 1.9
4	ZRH-CB-1816-3	Φ 3.0	10	1600	Φ 1.9
5	ZRH-CB-2416-3	Φ 3.0	10	1600	Φ 2.5
6	ZRH-CB-2418-4	Φ 4.0	10	1600	Φ 2.5
7	ZRH-CB-2423-3	Φ 3.0	10	2300	Φ 2.5
8	ZRH-CB-2423-4	Φ 4.0	10	2300	Φ 2.5

Disposable Ureteral Access Sheath

Features

- Hydrophilic coated guide sheath and dilator
- Metal reinforced wire to prevent guide sheath from bending and kinking
- Smooth joint, seamless integration, good insertion




Specification

No.	REF	Sheath I.D.		Working Length L ± 10(mm)
		Fr	mm	
1	ZRH-NQG-9.5-13	9.5	3.17	130
2	ZRH-NQG-9.5-20	9.5	3.17	200
3	ZRH-NQG-10-45	10	3.33	450
4	ZRH-NQG-10-55	10	3.33	550
5	ZRH-NQG-11-28	11	3.67	280
6	ZRH-NQG-11-35	11	3.67	350
7	ZRH-NQG-12-55	12	4.0	550
8	ZRH-NQG-13-45	13	4.33	450
9	ZRH-NQG-13-55	13	4.33	550
10	ZRH-NQG-14-13	14	4.67	130
11	ZRH-NQG-14-20	14	4.67	200
12	ZRH-NQG-16-13	16	5.33	130
13	ZRH-NQG-16-20	16	5.33	200

Disposable Nasobiliary Drainage Tube

Features

- Multiple shape selection
- Excellent bending resistance and shape memory performance
- Excellent X-ray development
- Multi side hole design, good drainage effect




Specification

REF	Stent Diameter		Applicable Channel	Working Length L ± 100(mm)
	Fr	mm		
A Type/Inverted A /Pigtail Type/Pigtail A/Pigtail/Inverted A /Short A	5	1.65	≥ 2.2	1700 (2600)
A Type/Inverted A /Pigtail Type/Pigtail A/Pigtail/Inverted A /Short A	6	1.98	≥ 2.8	1700 (2600)
A Type/Inverted A /Pigtail Type/Pigtail A/Pigtail	7	2.31	≥ 2.8	1700 (2600)
A Type/Inverted A /Pigtail Type/Pigtail A/Pigtail	8	2.64	≥ 3.2	1700 (2600)

Disposable Urinary Stone Retrieval Basket

Features

- Multiple specifications
- Unique handle design, easy to operate
- Headless end structure can be closer to the stone
- Multi-layer material outer tube
- 3 or 4 wire structure, easy to catch small stones




Specification

No.	Model	REF	Outer Sheath O.D. ± 0.1		Working Length L $\pm 10\%$ L (mm)	Opening Range(mm)
			Fr	mm		
1	3 wire	★ ZRH-WA-F1.712-8	1.7	0.56	1200	8
2		ZRH-WA-F1.712-15				15
3		★ ZRH-WA-F2.212-8	2.2	0.73	1200	8
4		ZRH-WA-F2.212-15				15
5		ZRH-WA-F312-8	3	1	1200	8
6		ZRH-WA-F312-15				15
7	4 wire	★ ZRH-WB-F1.712-10	1.7	0.56	1200	10
8		ZRH-WB-F1.712-15				15
9		★ ZRH-WB-F2.212-10	2.2	0.73	1200	10
10		ZRH-WB-F2.212-15				15
11		ZRH-WB-F312-10	3	1	1200	10
12		ZRH-WB-F312-10				15
13		ZRH-WB-F4.57-10	4.5	1.5	700	10
14		ZRH-WB-F4.57-15				15

Disposable Cleaning Brush

Features

- Disposable and reusable are available
- Unique brush design, easier to clean the endoscopic clamp hole and water vapor button channel
- Reusable cleaning brush: made of medical grade stainless steel, all metal, more durable
- Single head water and gas passage: used to clean water and gas button
- Double head water and gas passage: used to clean water and gas button
- Double brush head improves cleaning efficiency and cleanliness




Specification

Brush Type	REF	Channel Size Φ (mm)	Working Length L (mm)	Brush Diameter D (mm)
A Single Brush	ZHR-A-BR-0702	Φ 2.0	700 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
	ZHR-A-BR-1202	Φ 2.0	1200 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
	ZHR-A-BR-1602	Φ 2.0	1600 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
	ZHR-A-BR-2302	Φ 2.0	2300 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
B Double Brush	ZHR-B-BR-0702	Φ 2.0	700 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
	ZHR-B-BR-1202	Φ 2.0	1200 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
	ZHR-B-BR-1602	Φ 2.0	1600 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
	ZHR-B-BR-2302	Φ 2.0	2300 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
C Triple Brush	ZHR-C-BR-0702	Φ 2.0	700 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
	ZHR-C-BR-1202	Φ 2.0	1200 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
	ZHR-C-BR-1602	Φ 2.0	1600 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
	ZHR-C-BR-2302	Φ 2.0	2300 \pm 50	Φ 2.0/3.0/4.0/5.0/6.0
D Short Brush with Handle	ZHR-D-BR-0510	/	500 \pm 50	Φ 10.0

Disposable Polypectomy Hot Snare

Features

- Rope knitting technology, polyp is not easy to slip off
- Compatible with all mainstream high-frequency electric devices on the market
- Various specifications can be customized
- High-strength medical stainless steel, with fast cutting and good coagulation effect
- Handle and the snare can rotate synchronously to improve the operation efficiency and save time



Specification

No.	Snare Shape	Snare Shape D ^{-20%} (mm)	Working Length L ± 10% (mm)	Sheath O.D. d ± 0.1(mm)	Rotation Function (Yes / No)
1	Oval (A)	15/25	1200	∅ 1.8	No
2		15/25/35	1800	∅ 2.4	No
3		15/25/35	2300	∅ 2.4	No
4		15/25/35	1800	∅ 2.4	Yes
5		15/25/35	2300	∅ 2.4	Yes
6	Hexagonal (B)	15/25	1200	∅ 1.8	No
7		15/25/35	1800	∅ 2.4	No
8		15/25/35	2300	∅ 2.4	No
9		15/25/35	1800	∅ 2.4	Yes
10		15/25/35	2300	∅ 2.4	Yes
11	Semilunar (C)	15/25	1200	∅ 1.8	No
12		15/25	1800	∅ 2.4	No
13		15/25	2300	∅ 2.4	No

Disposable Hemostasis Clips

Features

- Short retention in the body, safer operation
- Large opening design of the clip, meet various treatment requirements
- It can be opened and closed repeatedly for repositioning
- 360 degree rotary design, convenient for clinical operation
- Stiffener design makes the clamp tighter, wide clip design makes the clamping more



Specification

No.	REF	Rotation Function YES / NO	Effective working length L(±100)mm	Clip max. Span W(+0.2/-2.7) mm	Instrument D(±0.15) mm
1	ZRH-HCA-235-15	YES / NO	2350	15	2.6
2	ZRH-HCA-195-15	YES / NO	1950	15	2.6
3	ZRH-HCA-165-15	★ YES / NO	1650	15	2.6
4	ZRH-HCA-235-12	★ YES / NO	2350	12	2.6
5	ZRH-HCA-195-12	YES / NO	1950	12	2.6
6	ZRH-HCA-165-12	YES / NO	1650	12	2.6

Disposable Injection Needle

Features

- Convenient single handed operation
- Various specifications apply widely
- Needle tip angle 30 degree, sharp puncture
- Blood can be observed by transparent inner tube



Specification

No.	REF	Sheath O.D. D ± 0.1 (mm)	Working Length L ± 50(mm)	Needle Gauge G (mm)	Needle Extension Length 1 ± 1 (mm)
1	ZRH-PN-2418-214	Φ 2.4	1800	21(0.800—0.830)	4
2	ZRH-PN-2418-234	★ Φ 2.4	1800	23(0.600—0.673)	4
3	ZRH-PN-2418-254	Φ 2.4	1800	25(0.500—0.530)	4
4	ZRH-PN-2418-216	Φ 2.4	1800	21(0.800—0.830)	6
5	ZRH-PN-2418-236	Φ 2.4	1800	23(0.600—0.673)	6
6	ZRH-PN-2418-256	Φ 2.4	1800	25(0.500—0.530)	6
7	ZRH-PN-2423-214	Φ 2.4	2300	21(0.800—0.830)	4
8	ZRH-PN-2423-234	★ Φ 2.4	2300	23(0.600—0.673)	4
9	ZRH-PN-2423-254	Φ 2.4	2300	25(0.500—0.530)	4
10	ZRH-PN-2423-216	Φ 2.4	2300	21(0.800—0.830)	6
11	ZRH-PN-2423-236	Φ 2.4	2300	23(0.600—0.673)	6
12	ZRH-PN-2423-256	Φ 2.4	2300	25(0.500—0.530)	6

Disposable Polypectomy Cold Snare

Features

- Rope knitting technology, polyp is not easy to slip off
- Various specifications can be customized
- High-strength medical stainless steel, with fast cutting and good coagulation effect
- Smooth and straight outer tube, convenient intubation



Specification

No.	Snare Shape	Snare Diameter D ^{-20%} (mm)	Working Length L ± 10% (mm)	Sheath O.D. d ± 0.1(mm)
1	Oval (A)	15/25	1200	Ø 1.8
2		15/25/35	1800	Ø 2.4
3		15/25/35	2300	Ø 2.4
4	Hexagonal (B)	15/25	1200	Ø 1.8
5		15/25/35	1800	Ø 2.4
6		15/25/35	2300	Ø 2.4
7	Semilunar (C)	15/25	1200	Ø 1.8
8		15/25	1800	Ø 2.4
9		15/25	2300	Ø 2.4